

NSF/ANSI STANDARD 49

Date Adopted & Effective	Issue Statements/ Annexes
<u>September 2026</u>	<u>199, 200, 203, 204, 205, 207, 208, 209, 211</u>
<u>November 2024</u>	<u>172, 173, 174, 176, 177, 179, 180, 181, 182, 186, 187, 190, 191, 192, 194, 195, 197, 198</u>
<u>March 2023</u>	<u>135, 142, 143, 161, 162, 166, 167, 168, 170, 171</u>
<u>March 2021</u>	<u>130, 141, 151, 153, 154, 155, 156, 157, 159, 160</u>
<u>March 2020</u>	<u>54, 82, 92, 120, 122, 127, 130, 133, 136, 138, 139, 140, 146, 147, 148, 149</u>
<u>January 2019</u>	<u>47, 59, 77, 105, 108, 109, 110, 111, 112, 115, 117, 118, 121, 125</u>
<u>March 2017</u>	<u>45, 56, 73, 76, 78, 79, 81, 86, 88, 90, 96, 99</u>
<u>February 2015</u>	<u>48, 49, 50, 51, 52, 53, 55, 60, 61, 72</u>
<u>July 2012</u>	<u>44, 46</u>
<u>November 29, 2010</u>	<u>23, 29</u>
<u>September 16, 2010</u>	<u>23, 24, 37/38, 41</u>
<u>June 2009</u>	<u>Annex G, 15, 28, 30, 34, 35, 36</u>
<u>April 2008</u>	<u>12, 13, 14, 16, 17, 18, 19, 20, 21, 22, 25, 26, 31</u>

Presented by James T. Wagner
jimwagner@cenvironment.com

Changes to NSF/ANSI 49

September 2026 changes

The following is based on the changes outlined in the Foreword of NSF/ANSI 49
The topics were discussed in committee as issue numbers which form the basis of discussion.

Issue 199, 200, 203, 204, 205, 207, 208, 209, 211

Only completed issues are discussed.

Issue 199

This revision adds language regarding the exhaust canopy connection setup instructions

- 5.22.4 Canopy connections listed as acceptable options for a BSC shall have a manufacturer specified set up instructions separate from the BSC's inflow and downflow setpoints, to ensure proper setup and function in the field.

SME's Take: The industry expectation is that proper setup is verified at each certification in addition to the canopy containment test.

Issue 199 (cont.)

This revision adds language regarding the exhaust canopy connection setup instructions

- N-1.13.2 indicates that manufacturers must establish setup instructions for their canopy.
- N-5.7.3.2.2.1
When possible, airflow balancing of the BSC being certified shall be completed following manufacturer instructions for canopy connection set up, including, but not limited to, bringing slot velocities or static duct pressure within the manufacturer specified ranges. The standard acknowledges this is not always possible. For example, installation of a canopy connection in a room with low ceiling clearance may make measurement of the exhaust static pressure impossible.

Issue 200

This revision updates and adds language regarding the effective cleaning and disinfection below the BSC workzone in section I-1.7.8

Use the following procedure to effectively clean or surface disinfect the area under the BSC work zone:

- 1) Don appropriate PPE.
- 2) Turn on blower and allow 5 min to purge the air.
- 3) Raise the sliding sash window to a full-open position.
- 4) Silence the audible alarm during the cleaning process.
- 5) After cleaning or surface disinfecting the BSC work zone surfaces, raise the front of the work surface so it is at an angle with the rear edge still supported by the cabinet. For larger cabinets, it may be necessary to have assistance or add a temporary support to hold the work surface at the angle needed for access.

Issue 200 (cont.)

This revision updates and adds language regarding the effective cleaning and disinfection below the BSC workzone in section I-1.7.8

- 6) Liberally spray or douse the area under the work surface with an appropriate surface decontaminant and let it sit for the required contact time for effectiveness.
- 7) Inspect the area for hazardous debris including broken glass, needles and other sharps. Take appropriate steps to safely remove this contaminated material.
- 8) After removal of hazardous debris and the contact time required for the decontaminant has elapsed, gently wipe the drain trough area. This gentle clean is intended to discover any hazardous debris that was not visible or was missed. After appropriate removal of any debris, continue to clean and disinfect the area under the work surface.

Issue 200 (cont.)

This revision updates and adds language regarding the effective cleaning and disinfection below the BSC workzone in section I-1.7.8

9) Inspect the drain valve. Over time debris and material may collect and block the drain valve. Inspect it to assure it is not blocked and will be functional if needed to manage and drain a spill.

10) Some BSCs have a paper catch screen installed to the rear of the area under the work surface in the return air path to the motor/blower. If the airflow is obstructed, performance of the BSC can be compromised. Therefore, the paper catch screen should be checked and obstructing material removed as appropriate. If paper products are used for procedures, the paper catch screen may require more frequent inspection. Removed paper must be properly discarded as contaminated hazardous waste.

Issue 203

This revision updates and adds language regarding the HEPA filter leak test in sections 5.19, 6.3, N-1.2.3.2, N-5.5.3.2, AND I-1.3.1.4

It also adds a normative reference for ANSI / ASSP Z9.14-2020 in Section 2.

ADDED LANGUAGE FOR TESTING HEPA FILTERS THAT CANNOT BE SCANNED

N-1.2.3.2 & N-5.5.3.2 (Added text)

When a cabinet is ducted so that the exhaust filter cannot be scanned, it may be leak tested by drilling a hole approximately 0.3 in. (8 mm) in diameter in the duct at a downstream location that will produce a well-mixed aerosol, and inserting the photometer sampling probe with rigid extension tubing through the hole **approximately eight duct diameters 8 or more from the HEPA filter or according to manufacturer recommendations.**

SME's Take: Leaks may be missed if probe testing too close to the HEPA filter in question

Issue 203 (cont.)

This revision updates and adds language regarding the HEPA filter leak test in sections 5.19, 6.3, N-1.2.3.2, N-5.5.3.2, AND I-1.3.1.4

It also adds a normative reference for ANSI / ASSP Z9.14-2020 in Section 2.

- **5.19** ...If filters cannot be scan tested, they shall be tested for total leakage not to exceed 0.005% in accordance with Section N-1.2
- **6.3.3** When filters cannot be scan tested, polydisperse DOP or equivalent sustained penetration shall not exceed 0.005% of the upstream concentration at any point when measured on a linear or logarithmic scale photometer in a downstream location that will produce a well-mixed aerosol approximately eight duct diameters or more from the HEPA filter, or according to manufacturer recommendations (see high-efficiency air filters).

Issue 203 (cont.)

This revision updates and adds language regarding the HEPA filter leak test in sections 5.19, 6.3, N-1.2.3.2, N-5.5.3.2, AND I-1.3.1.4

It also adds a normative reference for ANSI/ASSP Z9.14-2020 in Section 2.

I-1.3.1.4 (Added text)

- Type B BSCs require a higher static pressure that must increase as their exhaust filters load. They must be direct hard connected to the on an exhaust duct and fan dedicated to that individual cabinet, and not be ganged with other Type B BSCs, or other exhausted laboratory devices requiring a lower static pressure (e.g. fume hoods, canopy-connected BSCs) **and should have duct access approximately eight duct diameters 8 from the HEPA filter or more, or as per manufacturer recommendations downstream of the cabinet for exhaust filter integrity testing**

Issue 204

This revision updates and adds language regarding cabinet types in sections N-5.3.3.2.2 and N-5.3.3.2.4

N-5.3.3.2.2 for constricted access inflow velocity reading using a thermal anemometer now applies to C1 and B1 cabinets as well as A1, A2, and B2 cabinets.

Issue 204 (cont.)

This revision updates and adds language regarding cabinet types in sections N-5.3.3.2.2 and N-5.3.3.2.4

N-5.3.3.2.4

Canopy-connected A1, A2, and C1 cabinets shall be tested with a method that measures the inflow volume at the work access opening.

Issue 205

This revision updates and adds language regarding the alternative inflow measurement methods in sections N-5.3.3.2

Minimum Clearance to use a DIM

2024...the space between the face of the BSC and permanent fixture directly opposite the access opening is less than 42 in (1.1 m);

2026...the space between the base (intake) of the DIM and permanent fixture directly opposite the access opening is less than 6 in. (150 mm)

Issue 207

This revision updates *Bacillus subtilis* to *Bacillus atrophaeus* throughout the standard.

Issue 208

This revision updates and adds language regarding the light intensity test in sections N-1.4.3 and N-5.9.3 to match language used in Figure 13.

N-1.4.3 & N-5.9.3

Measure the cabinet Measure the cabinet light intensity along the ~~side-to-side~~ **front-to-rear** centerline of the work tray on a uniform linear pattern in increments close to but not greater than 12 in (300 mm), starting 6 in (150 mm) from the side walls.

SME's Take: Actual procedures don't change; language was adapted to fit what is actually done.

Issue 209

This revision updates and adds language regarding the pressure testing in sections N-1.1.1.2 and N-1.11.2

N-1.1.1.2 & N-1.11.2 (Added text)

- manometer, pressure gauge, or pressure transducer system with a minimum range of 0 to 2 in. w.g. (0 to 500 Pa) and accurate to $\pm 2\%$ of reading ± 0.001 in. w.g. (0.2 Pa) **or $\pm 1\%$ of reading ± 0.01 in. w.g. (2 Pa);**

Issue 211

This revision updates and adds language regarding the B1 and B2 exhaust alarm language in section N-5.7.3.2.1

N-5.7.3.2.1

reduce exhaust volume by ~~at least~~ **no more than** 20% once the cabinet is set or certified in its acceptable air-flow range.

SME's Take: Actual procedures don't change; language was adapted to fit what is actually done.

Issue 211 (cont.)

This revision updates and adds language regarding the B1 and B2 exhaust alarm language in section N-5.7.3.2.1

~~NOTE — For direct connected Type B1 or B2 BSCs, measure the static pressure in the duct-work between the cabinet and duct-mounted balancing dampers.~~

- measured static pressure in the ductwork between the cabinet and duct-mounted balancing damper, if applicable
- measured static pressure in the ductwork above the cabinet when there is no balancing damper, if applicable.

SME's Take: A measurement of duct static pressure is now considered part of mandatory tests and no longer a note of recommendation.

Changes to NSF/ANSI 49

November 2024 changes

The following is based on the changes outlined in the Foreword of NSF/ANSI 49
The topics were discussed in committee as issue numbers which form the basis of discussion.

Issue 172, 173, 174, 176, 177, 179, 180, 181, 186, 187, 190, 191, 192, 194, 195, 197, 198

Only completed issues are discussed.

Issue 172

This revision updates language regarding the pressure decay test in Section 6.2 and Sections N-1.1.1.3 through N-1.1.1.5.

- 6.2 Pressure decay / soap bubble / tracer gas leak
- 6.2.1
 - ~~The cabinet shall hold 2 in w.g. (500 Pa) within +/- 10% for 10 min.~~
 - The cabinet shall hold 0.5 in w.g. (125 Pa) for 10 min with decay not exceeding 10%.
- 6.2.2
 - For manufacturer testing only, the soap bubble method may be used when pressure plates fail: all welds, gaskets, penetrations, or seals on exterior surfaces of air plenums shall be free of soap bubbles when at ~~2 in w.g. (500 Pa)~~ 0.5 in w.g. (125 Pa) $\pm$ 10% pressure above atmospheric.
- N-1.1.1.3 Method (pressure decay)
 - d) Pressurize the cabinet with air to a reading of ~~2 in w.g. (500 Pa)~~ 0.5 in w.g. (125 Pa), turn off the pressurizing air, and measure the pressure after 10 min.
 - e) If the cabinet ~~does not hold pressure withing 10% after 10 min~~ loses > 10% of the initial pressure after 10 min, use the liquid leak detector to check for leaks in the pressure plates used to seal the access opening, exhaust, and fan inlet (where applicable)...

Issue 172 Cont.

- N-1.1.1.4 Method (soap bubble)
 - d) Pressurize the cabinet with air to ensure a continuous reading of ~~2 (500 Pa)~~ 0.5 (125 Pa) in w.g. $\pm 10\%$.
- N-1.1.1.5 Acceptance
 - N-1.1.1.5.1 Pressure decay
 - ~~The cabinet shall hold 2 in w.g. (500 Pa) $\pm 10\%$ for 10 min.~~ The cabinet shall not lose $> 10\%$ of the initial pressure after 10 min, when starting at 0.5 in w.g. (125 Pa). This requirement shall be met for all cabinet design certification testing.
 - N-1.1.1.5.2 Soap bubble
 - All welds, gaskets, penetrations, and seals on exterior surfaces of air plenums shall be free of soap bubbles when at ~~2 in w.g. (500 Pa)~~ 0.5 in w.g. (125 Pa) $\pm 10\%$ pressure above atmospheric. This requirement may be met during manufacturing as an alternative to the pressure decay test.

Issue 173

This revision adds definitions for **cleanable**, **easily cleanable**, and **tubing restraint** and affirms revised and new language related to design and construction requirements throughout Section 5. It also removes former Sections 5.7.1.3, 5.11, 5.12, 5.18.

- **3.14 cleanable:** Able to be freed of residues of soiling materials.
- **3.20 easily cleanable:** Manufactured so that soiling material may be removed by manual cleaning methods.
- **3.37 tubing restraint:** Device to minimize movement (i.e., $\pm 1/4$ in) of tubing.
- ~~5.7.1.3 **Fillet material** Parent material or hard solder may be used as fillet material in structurally sound seams~~
- ~~5.11 **Welds** Welds shall meet the smoothness requirements of the applicable surface.~~
- ~~5.12 **Solder** Solder shall only be used to seal structurally sound seams or as fillet material (see Section 5.7.1.3)~~
- ~~5.18 **Fixed Panels** Fixed panels shall be designed, constructed, and fastened to eliminate projections and openings.~~

Issue 173 Cont.

- 5.2 ~~Cleanability~~ Surface decontamination

- Interior work, exposed interior, and the other interior surfaces subject to splash or spillage shall be readily accessible and easily cleanable as assembled or when removed. ~~Interior work, exposed interior, and other interior surfaces, including plenums, shall be capable of being vapor or gas decontaminated.~~

- 5.3 Gas / vapor decontamination²¹

- ~~Cabinets shall be designed to be decontaminated~~ Interior work, exposed interior, and other interior surfaces, including plenums, shall be capable of being vapor or gas decontaminated with an inactivating agent (such as formaldehyde gas) without being moved. Closure to contain decontaminating agents should be limited to gas-tight sealing of air intake and exhaust openings with metal plates, or plastic film and tape, or equivalent.²¹
- ~~Pressure tight valves, if provided, suitable for decontamination shall be located on the clean side of the HEPNULPA filter.~~

- 5.75.9 Joints and seams

- 5.7.15.9.1 Interior work and exposed interior surfaces Total work area and exposed drain pan surfaces
...All other seams shall be sealed. Equipment parts shall be stamped, extruded, formed, or cast in one piece. Joints shall be fabricated to eliminate dirt-catching horizontal ledges.
- 5.7.25.9.2 Other interior and exterior surfaces
...All joints and seams subject to exposure to vapor, or toxic volatile substances, or both, and exposed to the outside environment shall be sealed. All other seams shall be closed.

Issue 173 Cont.

- 5.85.10 Fastening methods

- 5.8.15.10.1 Exposed fastenings

Exposed screw threads, projecting screws, and studs shall not be used on interior work surfaces. They shall only be used on exposed interior and other interior surfaces when other fastening methods are impractical. All metal fasteners and studs subject to maintenance shall not be subject to excessive overspray.

- 5.10.2 Exterior fastenings

Fasteners for exterior removable panels that are gasketed and subject to pressure shall be studs with solid acorn nuts, or equivalent, so that the gasket is sealed. Fasteners for other removable panels may be low profile type fasteners (truss, round counter sunk, flat counter sunk head [see Figure 4]), or studs with solid acorn nuts. All metal fasteners and studs subject to maintenance shall not be subject to excessive overspray.

- 5.10.32 Interior fastenings

In areas subject to cleaning, interior fastenings and joinings shall be fabricated to minimize projections, ledges, and recesses. All metal fasteners and studs subject to maintenance shall not be subject to excessive overspray.

Issue 173 Cont.

- 5.11 Removeable panels
 - ...The design and construction of removable panels shall minimize projections and openings. Removable panels for access into contaminated areas shall be designed so that, upon reassembly, a seal is provided without application of additional sealant as required in Section 6.2.5.20.
- 5.12 Stability
 - ...in accordance with Section N-1.7.
- 5.14 Legs and feet
 - Legs and feet shall be sufficiently rigid to provide support with a minimum of cross bracing. They shall be fastened to the cabinet and shaped at floor or bench top contact to minimize the accumulation of splash and spillage. Legs and feet shall be of simple design, with no exposed threads. ...

Issue 173 Cont.

- 5.20 Gaskets and sealants

...

— ~~HEPA/ULPA filter seals shall be leakproof when tested in accordance with Section N-1.3.~~
Gaskets on HEPNULPAHEPA/ULPA filters shall have Interlocking corners or sealed joints;

— ...

— gaskets used in cabinet seams or on the facing of service panels shall provide for resealing after removal and reinstallation;

— service panels providing access to contaminated spaces shall have a gasket; and

— ...

- 5.22.2 Internal cabinet supply/ exhaust fan interlock alarm

- ... Audible and visual alarms shall activate within 15 s of the alarm condition initiating.

- 5.22.6 Types A1, A2, or C1 cabinet low inflow alarm

... If the manufacturer-specified inflow velocity alarm set point is more than 10 ft/min (0.051 m/s) less than the nominal inflow velocity, the test as specified in Section N-1.6.3.1.h will be performed with the inflow velocity at this set point ± 3.0 (0.015 m/s) ± 2.0 ft/min (0.01 m/s).

~~If the manufacturer-specified inflow velocity alarm set point is no more than 10 ft/min (0.051 m/s) less than the nominal inflow velocity, the inflow alarm point shall be tested as specified in Section N-1.6.3.1.h.~~

Issue 173 Cont.

- 5.23 Electrical components
 - 5.23.1 Motor

- a thermal protector shall be provided. It shall not trip at 115% of the rated voltage under maximum load and ambient temperature conditions. The motor shall be rated for continuous operation;
- fan motors shall be sized to operate at a static pressure sufficient to meet the requirements of Section 6.13;
- All fan motors shall be variable speed and shall have controls that can be secured. Controls shall...
- motors and lights shall be separately protected from the receptacles. Circuit overload protection conforming to the National Electrical Code¹⁴ shall be provided. Flexible power cords for single-phase power shall be 3-wire, with the ground wire connected to the frame, unless otherwise specified and sized in accordance with the National Electrical Code¹⁴ for the specified load(s).

Issue 173 Cont.

- 5.25 5 Gauges, pressure transducers, and switches
 - Pressure gauges indicating the differential pressure across the recirculated air filter, if may be provided, but are not required. shall be installed in accordance with the manufacturer's instructions. Hose connections to the gauge and sampling port. All tubing shall have an inside diameter less than the outside diameter of the mating connection and shall be secured by one of the following methods when the tubing is coming from a potentially contaminated plenum:
 - positive compression clamps;
 - tubing restraint, securing the tubing within 4 in (10 cm) of each connection point;
 - a distance not > 1 in (2.5 cm) between connections points opposite each other, with the length of tubing exceeding the distance between the points and each connection point secured to the same solid surface.
 - ... Tubing coming from positively pressurized and potentially contaminated plenums shall be protected by an in-line HEPA filter, with a minimum rated efficiency of 99.97%.

Issue 173 Cont.

- 5.28 Total work area components placement
 - ... work surfaces, intake air grills, sash wipes, and exhaust air grills...
- 5.29.2 The following shall be readily viewable from the front of the cabinet
 - -software or firmware version number(s) used by the BSC.
- 5.31 User-modified pass-throughs
 - Retaining element(s) for the various cables, tubes, etc. shall be readily replaceable by the user;
 - The pass-through shall bear cautionary labels both interior and exterior referencing use; and
 - The pass-through shall be connected to a negative pressure zone to prevent contaminated air from escaping the cabinet.

Issue 174

This revision affirms revised language related to recommended preventative maintenance procedures in Section N-5.1.

- N-5.1 Field certification preconditions and intervals
 - Field certification may not address other recommended maintenance procedures for specific cabinets, stands, and accessories. Manufacturer recommended preventative maintenance procedures shall be provided in the cabinet's operating manual.

Issue 176

This revision affirms revised language related to the verification of the Secondary Inflow Method throughout Section N-1.8.3 and N-5.3.3.

Issue 176 Cont.

- N-1.8.3.1 General
 - The nominal set point average inflow velocity shall be determined by a direct inflow reading instrument measurement. After the nominal set point is determined by a direct inflow reading instrument measurement, readings shall be taken by the appropriate alternate calculated or measured method recommended by the manufacturer. Both of these set point values shall meet the requirements of Section N-1.9.4.
- N-1.8.3.4 Alternate inflow measurement methods
 - These methods, approved by the testing organization, shall be validated and provided by the manufacturer and shall be subject to review by the testing organization. Manufacturer validation procedures shall contain no fewer than ten replicate tests. The testing organization's approval shall be based on review of data and successful reproduction of test results. The following methods have been found to be acceptable on some cabinets:
- N-1.8.3.4.2.b
 - ...The air measurement probe shall be held rigidly in a freestanding fixture provided by the manufacturer that permits accurate positioning and does not distort the airflow pattern. This includes a ring stand and clamp or manufacturer supplied probe holder, or if specified by the BSC manufacturer, taping the probe to the inside or outside of the sash. The anemometer probe shall not be hand held.
- N-1.8.3.4.5 Alternate method acceptance
 - Acceptance criteria of the alternate method shall be based on inflow determined by the direct measurement. The fully corrected alternate method inflow velocity shall be within ± 5 ft/min (± 0.025 m/s) of the direct measurement results when both measurements are completed on the same day.

Issue 176 Cont.

- **N-5.3.3 Methods**

- One of these methods was validated per cabinet model and provided by the manufacturer, which was reviewed and approved by the testing organization. ~~Manufacturer validation procedures contained no fewer than ten replicate tests. The testing organization's approval will be based on review of data and successful reproduction of test results. The validated alternate method is on the manufacturer's data plate.~~

- **N-5.3.3.1 General**

- ~~When the testing organization has determined the nominal set point on a given model and size of cabinet using a DIM device, and an appropriate alternative method has been validated for that cabinet by the testing organization, this alternate method may be used to establish the set point on the same model and size of cabinet in the field.~~

Issue 177

This revision affirms revised and new language related to noise level tests in Sections 6.4.2, N-1.3.3.c, N-1.3.4, N-5.11.3.c, N-5.11.3.g, and N-5.11.4.

Issue 177 Cont.

- 6.4 Noise level
- 6.4.2
 - The overall noise level ~~12 in (300 mm)~~ 14 in (360 mm) in front of the top of the access opening and ~~15 in (380 mm)~~ above the plane of the work surface at the vertical centerline and 4 in (100 mm) above the top of the access opening at the side-to-side centerline of the cabinet shall not exceed:
 - 67 dBA with a maximum background level of 57 dBA (biosafety cabinets with an access opening less than 12 in [300 mm]); and
 - 69 dBA with a maximum background level of 59 dBA (biosafety cabinets with an access opening greater or equal to 12 in [300 mm]).

Issue 177 Cont.

- N-1.3 Noise level test
- N-1.3.3 Method
 - c) Position the noise level meter ~~12 in (300 mm)~~ 14 in (360 mm) in front of the sash at the top of the access opening and ~~15 in (380 mm)~~ above the plane of the work surface 4 in (100 mm) above the top of the access opening, in line with the vertical side-to-side centerline of the cabinet (see Figure 12).
- N-1.3.4 Acceptance
 - Biosafety cabinets with an access opening:
 - < 12 in (300 mm): The net noise level in front of the cabinet shall not exceed 67 dBA; or
 - ≥ 12 in (300 mm): The net noise level in front of the cabinet shall not exceed 69 dBA.

Issue 177 Cont.

- N-5.11 Noise level tests
- N-5.11.3 Method
 - c) Position the noise level meter ~~12 in (300 mm)~~ 14 in (360 mm) in front of the ~~cabinet front edge~~ sash at the top of the access opening and ~~15 in (380 mm)~~ above the plane of the work surface 4 in (100 mm) above the top of the access opening, in line with the vertical side-to-side centerline of the cabinet (see Figure 12).
 - ...
 - g) Reported values shall be:
 - gross sound level reading (cabinet sound level before correcting for background);
 - background sound level reading (with cabinet turned off);

Issue 177 Cont.

- N-5.11.4 Acceptance
 - Overall noise level in front of the cabinet shall not exceed 70 dbA (72 dbA for 12 in opening cabinets) when measured where the maximum ambient sound level is no > 60 dbA. When the ambient sound level is > 60 dbA, the reading obtained in Section N-5.11.3.c shall be corrected in accordance with curves or tables provided in the instrument operator's manual. If this information is not available, standard correction curves or tables shall be used (see Table N-5.1).

Biosafety cabinets with an access opening < 12 in (300 mm): The ~~net~~ cabinet noise level in front of the cabinet shall not exceed 70 dBA.

Biosafety cabinets with an access opening $\geq$ 12 in (300 mm): The cabinet noise level in front of the cabinet shall not exceed 72 dBA.

Table N-5.1
Correction chart for sound level readings

Difference between total and background sound readings in dbA	Number to subtract from total to yield corrected noise level
0 to 2	reduce background levels
3	3
4 to 5	2
6 to 10	1
> 10	0

Issue 179

- This revision affirms new language related to the requirements for canopy connections as new Section 5.29.3.
- 5.29 Data plate(s)
- 5.29.3
 - A canopy connection is not required. However, when a canopy is submitted for approval as an acceptable option, it shall have a physical data plate including all of the following information: — manufacturer's name and address; and — canopy model number.

Issue 180

This revision corrects language regarding motor blower performance procedure Section N-1.11.3.c.

- N-1.11 Moter / blower performance
- N-1.11.3 Method
 - c) ... The manufacturer shall locate the negative pressure tap not less more than one-half equivalent diameter from the blower inlet...

Issue 181

This revision updates and adds language regarding motor blower performance procedure Sections N-1.11.2 and N-1.11.3.b.

- N-1.11.2 Apparatus
 - Updated references
- N-1.11.3 Method
 - b) Measure the total airflow volume rate, ft^3/min (m^3/s), and determine that the cabinet blow is delivering at the nominal set point (see Sections N-1.8 and N-1.9). The cabinet supply air volume shall be determined as in section N-1.9.3.4.4.d. The total airflow volume rate is the flow from all cabinet fans and does not include facility exhaust. The downflow volume rate is always measured for this test. The inflow volume rate is measured only for Type A and Type C cabinets. Measure inflow volume rate following the procedure in Section N-1.8, using the direct inflow reading instrument. Measure the downflow velocity following the procedure in Section N-1.7 and then multiply the overall average downflow velocity by the area of the measurement plane to calculate the downflow volume rate. For Type A and Type C cabinets, add the inflow and downflow volume measurements together for total airflow volume.

Issue 182

- This revision updates language related to the requirements canopy connection acceptance criteria in Section N-1.13.5.
- N-1.13 Canopy connection test
- N-1.13.5
 - Inflow velocity shall not be reduced by more than 10 ft/min (0.051 m/s) to less than the lowest value verified by the NSF/ANSI 49 biological challenge testing after turning the facility exhaust off.

Issue 186

- This revision removes former Sections 6.17.1 and N-1.13.3.1 related to drop testing.
- **6.17 Air velocity stability**
- ~~6.17.1~~
 - ~~When the cabinet is subjected to a 10-mm free fall drop on each side, the cabinet inflow velocity and downflow velocity (where applicable) shall not change by more than 5 ft/min (0.025 m/s). There shall be no visible damage to the cabinet following the shock.~~
- ~~N-1.13.3.1 Shock Stability~~

Issue 187

This revision updates power failure disconnection time from 1 h to 5 min in Section 6.17.2 and Section N-1.12.3.2.

- 6.17 Air velocity stability
- 6.17.2
 - When the cabinet has been disconnected from power for a minimum of ~~1 hour~~ 5 min, the cabinet inflow velocity (where applicable), or downflow velocity (where applicable), or both, shall not change by more than 3 ft/min (0.015 m/s) when power is restored...
- N-1.12.3.2 Power failure stability
 - c) Disconnect power to the cabinet for a minimum of ~~1 hr~~ 5 min.

Issue 190

This revision updates language related to the personal protection test in Section N-1.6.3.1.

- N-1.6.3 Personnel protection test
- N-1.6.3.1 Method
 - f) Filter the sampling fluid from all of the AGI-30 samplers³⁵ through a 47-mm diameter ~~0.2- to 0.22-µm~~ 0.45-µm membrane filter, remove the filter aseptically, and place it on appropriate media...

Issue 191

This revision affirms revised language related to control plate placement in Sections N-1.6.3.1.c and N-1.6.4.1.e.

- N-1.6.3 Personnel protection test
- N-1.6.3.1 Method
 - c) ... As a positive control, an agar plate shall be placed under the center of the cylinder, and supported a minimum of 0.50 in (13 mm) above work surface or below the front intake grill, to minimize in a manner that minimizes the obstruction of airflow into the grill (see Figures 16 and 17).
- N-1.6.4 Product protection test
- N-1.6.4.1 Method
 - e) As a positive control, an agar plate shall be placed under the center of the cylinder, and supported a minimum of 0.50 in (13 mm) above or below the front intake grill, to minimize the work surface or front intake grill in a manner that minimizes the obstruction of airflow into the grill (see Figure 19).

Issue 192

This revision updates a reference in Section 6.8.1 and removes previous Sections N-1.7 through N-1.7.2 relating to the resistance to overturning test.

- 6.8.1 Resistance to overturning
 - Cabinets shall conform to the requirements of ~~UL 61010-1¹⁷ or current edition,~~ Section 7.3 IEC 61010-1.¹⁷
- ~~N-1.7 Resistance to overturning~~
 - ~~N-1.7.1 Method~~
 - a) ~~Block the cabinet (adjusted to manufacturer's tallest rated service position on stand if applicable) at front or rear bottom edge to prevent lateral movement.~~
 - b) ~~Tilt the cabinet 10° from horizontal in the direction most likely to cause overturn.~~
 - ~~N-1.7.2 Acceptance~~

~~The cabinet shall not initiate overturn when tilted 10° from horizontal in the direction most likely to cause overturn~~

Issue 194

This revision updates language related to the noise level acceptance criteria in Section N-5.11.4 and adds Table N-5.1.

Issue 194 Cont.

- N-5.11.4 Acceptance

- Overall noise level in front of the cabinet shall not exceed 70 dbA (72 dbA for 12 in opening cabinets) when measured where the maximum ambient sound level is no > 60 dbA. When the ambient sound level is > 60 dbA, the reading obtained in Section N-5.11.3.c shall be corrected in accordance with curves or tables provided in the instrument operator's manual. If this information is not available, standard correction curves or tables shall be used (see Table N-5.1).

Biosafety cabinets with an access opening < 12 in (300 mm): The net cabinet noise level in front of the cabinet shall not exceed 70 dBA.

Biosafety cabinets with an access opening $\geq$ 12 in (300 mm): The cabinet noise level in front of the cabinet shall not exceed 72 dBA.

Table N-5.1
Correction chart for sound level readings

Difference between total and background sound readings in dbA	Number to subtract from total to yield corrected noise level
0 to 2	reduce background levels
3	3
4 to 5	2
6 to 10	1
> 10	0

Issue 195

This revision updates language related to the installation and lifespan of BSCs throughout Sections I-1.3, I-1.8, and I-1.9.

Issue 195 Cont.

- I-1.3.3.5 Base stands, supports, restraints
 - Base stands with and without lifts or supports or required restraints (i.e., seismic) shall be considered at the time of specification. Some models of cabinets can weigh up to 900 lb (408 kg). The BSC must be attached to a manufacturer recommended base stand to support the unit's weight. All base stands have a maximum height specified by the manufacturer to prevent overturning of the BSC; this maximum should never be exceeded. Seismic restraints, if required, shall meet all state and local building codes and receive approval from the local authority having jurisdiction.

Issue 195 Cont.

- I-1.3.3.6 Mobile or non-permanent installations
 - BSCs used with mobile or castored base stands with and without lift capability have been used when the BSC is lifts can be operated in multiple routine marked locations in the same or adjoining laboratories. Proper cabinet operation airflow performance should be confirmed by airflow pattern testing at each site of use marked routine location used. If the cabinet BSC on a mobile or castored base stand is relocated to another non-routine marked location, laboratory, or facility, or and is subjected to excessive shock, or vibration, or both, during the move, the BSC should be recertified and airflow performance again checked at each location used to ensure it is functioning in a proper manner.

Issue 195 Cont.

- I-1.8 Moving a permanently installed BSC
- I-1.8.2
 - ... Depending on circumstances of the move, (i.e., access for service around cabinet, cabinet use, new permanent location, etc.), BSCs may be required to be space decontaminated before the move. After a BSC is moved, it ~~should~~ may also be required to be re-certified according to applicable performance standards.
- I-1.9 BSC lifespan
 - ... BSCs manufactured in the 1970s, 1980s, 1990s and through the early 2000s have provided over 15 yr of service...
 - ... T12 fluorescent bulbs may not be available, as well as single phase AC motors...

Issue 197

This revision updates and adds new language related to chemical resistance testing in Sections 4.5.2.1, N-4.1.1, and N-4.1.2.

- 4.5.2.1 Chemical resistance
 - ... ~~5~~2% phenol; and
- N-4.1.1 Chemicals
 - The following chemicals shall be used for resistance testing:
 - 1N hydrochloric acid, contact time = 4 h;
 - 1N sodium hydroxide, contact time = 4 h;
 - 1% quaternary ammonium compound; , contact time = 4 h;
 - 5% formaldehyde, contact time = 4 h;
 - 5,000 ppm hypochlorite, contact time = 4 h;
 - 2% iodophor, contact time = 4 h;
 - ~~5~~2% phenol, contact time = 2 h; and
 - 70% ethyl alcohol (ethanol), contact time = 4 h.
- N-4.1.2 Method
 - Reagents shall be allowed to remain on the surface for ~~4 hr~~ contact time specified in Section N-4.1.1, and tests shall be performed so the testing surface is wet throughout the entire test period. After ~~4 hr~~ contact time specified in Section N-4.1.1, the surface shall withstand...

Issue 198

This revision adds new language related to replacement filters in Section 5.19.

- 5.19 Filters

...

- HEPA/ULPA filters shall be specified (and considered equivalent) as to:
- type or grade (i.e., per IEST-RP-CC00112);
- external physical dimensions (width × depth × height);
- location and description, including the material of the gasketing on the filter;
- presence and description of protective screens;
- other structures on or inside the filter, specifically structural elements and media separator type;
- material and construction of filter frame;
- internal physical dimensions of the media pack unobstructed by potting or sealing media; and
- pressure drop across the filter ± 0.05 in w.g. (12.5 Pa) at the rated flow.

...

CHANGES TO NSF/ANSI 49: MARCH 2023 CHANGES

Changes Effective March 2023

The following is based on the changes outlined in the Foreword of NSF/ANSI 49

The topics were discussed in committee as issue numbers which form the basis of discussion.

- Issue 135, 142, 143, 161, 162, 166, 167, 168, 170, 171
- Only completed issues are discussed.

ISSUE 135

This revision added language regarding the power failure stability for Type B and C1 cabinets in Section N-1.13.4.3.

N-1.13.4.3 Power failure stability

The difference between the initial inflow velocity and the final inflow velocity shall not exceed 3 ft/min (0.015 m/s). The difference between the initial downflow velocity and the final downflow velocity shall not exceed 3 ft/min (0.015 m/s). The cabinet lights shall come back on automatically when power is restored. The cabinet blower shall come back on automatically when power is restored, **except in the cases of all Type B or those C1 cabinets connected to an exhaust system in which sensors determine when there is insufficient exhaust flow.** Alarm parameters (if so equipped) shall remain unchanged from the set points prior to power loss. The cabinet shall provide the user with a visual indication that there was a loss of power.

*This section was **modified**. The **highlighted** language was added to the Power failure stability requirement to prevent cabinets shut down because of insufficient exhaust from automatically restarting.*

ISSUE 142

This revision adds language regarding field tests related to specific industry use of BSCs as Section N-5.1.3.

N-5.1.3 Tests related to specific industry use

When field testing any Class II BSC, the certifier should take into consideration any additional testing required by a specific industry use. Examples include but are not limited to pharmaceuticals and electronics. Any additional testing completed is an adjunct to, and not a substitute for, the minimum testing required to consider a BSC field tested in accordance with NSF/ANSI 49.

*This section **was added** to accommodate differences in applications of BSC from industry to industry. For example, the sterile compounding industry specifically requires particle counting to be done in all Primary Engineering Controls (PECs) which includes Class II BSCs. When certifying to USP Chapter <797> tests beyond those needed for compliance with NSF/ANSI standard 49 are required but those tests must be done in addition to the NSF tests, not in place of.*

ISSUE 143

This revision modifies language regarding performance testing of Noise Level in Sections N-1.3 and N-5.11.

N-1.3 Noise level test

N-1.3.1 Purpose

This test provides a uniform method for measuring the noise level produced by the cabinet. The methods can be performed in most acoustically ordinary rooms, such as a factory, where walls are neither sound absorbing nor completely sound reflecting. The cabinet shall be operated at the manufacturer's recommended nominal set points **± 2 ft/min (0.01 m/s).**

This section was modified from +/- 3 ft/min to +/- 2 ft/min.

ISSUE 143 (CONTINUED)

This revision modifies language regarding performance testing of Noise Level in Sections N-1.3 and N-5.11.

N-1.3 Noise level test

N-1.3.1 Purpose

This test provides a uniform method for measuring the noise level produced by the cabinet. The methods can be performed in most acoustically ordinary rooms, such as a factory, where walls are neither sound absorbing nor completely sound reflecting. The cabinet shall be operated at the manufacturer's recommended nominal set points **± 2 ft/min (0.01 m/s).**

This section was modified from +/- 3 ft/min to +/- 2 ft/min.

ISSUE 143 (CONTINUED)

This revision modifies language regarding performance testing of Noise Level in Sections N-1.3 and N-5.11.

N-1.3.3 Method

- a) Turn on the cabinet blower and lights.
- b) Set the instrument to the "A" weighting mode.
- c) **Position the noise level meter** 12 in (0.30 m) in front of the cabinet **front** edge of the access opening and 15 in (0.38 m) above the plane of the work surface, in line with the vertical centerline of the cabinet (see Figure 12).
- d) Measure the gross noise level.**
- e) Measure the background noise level with the cabinet blower(s) and light(s) off and, if applicable, the exhaust blower on.**
- f) Correct the gross noise level in accordance with curves or tables provided in the instrument operator's manual to determine the net noise level.**

N-1.3.4 Acceptance

The net noise level in front of the cabinet shall not exceed 67 dBA.

Item c) was modified for clarity.

Items d),e),f) and N-1.3.4 Acceptance were rewritten or added for clarity

Note that the process and acceptance criteria did not change.

ISSUE 143 (CONTINUED)

This revision modifies language regarding performance testing of Noise Level in Sections N-1.3 and N-5.11.

N-5.11.3 Method

- a) Operate the cabinet within 5 ft/min (0.025 m/s) of the nominal set point with lights on.
- b) Set the instrument to the "A" weighting mode.
- c) Position the noise level meter 12 in (300 mm) in front of the cabinet (front edge of the access opening) and 15 in (380 mm) above the plane of the work surface, in line with the vertical centerline of the cabinet (Annex N-1, Figure 13).

d) Measure the gross noise level.

e) Measure the background noise level with the cabinet blower(s) and light(s) off and, if applicable, the exhaust blower on.

f) Correct the gross noise level in accordance with curves or tables provided in the instrument operator's manual to determine the net noise level.

g) Reported values shall be:

- gross sound level reading;
- background sound level reading;
- net sound level;
- pass or fail; and
- name of test (noise level tests).

*Highlighted items d),e),f) and g) were rewritten for clarity
Note that the process did not change.*

ISSUE 143 (CONTINUED)

This revision modifies language regarding performance testing of Noise Level in Sections N-1.3 and N-5.11.

N-5.11.4 Acceptance

The net noise level in front of the cabinet shall not exceed 70 dBA.

The new language above replaced the old language below. The language was changed for clarity, but the effective acceptance criteria did not change.

N-5.11.4 Acceptance

Overall noise level in front of the cabinet shall not exceed 70 dbA when measured where the maximum ambient sound level is no greater than 60 dbA. When the ambient sound level is greater than 60 dbA, the reading obtained in Section N-5.11.3.c) shall be corrected in accordance with curves or tables provided in the instrument operator's manual. If this information is not available, standard correction curves or tables shall be used (see following table).

ISSUE 143 (CONTINUED)

This revision modifies language regarding performance testing of Noise Level in Sections N-1.3 and N-5.11.

N-5.11.4 Acceptance

The noise correction table below was eliminated from this version of NSF/ANSI 49 in favor of correction curves found in the instrument operators manual.

Difference between total and background sound readings in dbA	Number to subtract from total to yield corrected noise level
0 to 2	reduce background levels
3	3
4 to 5	2
6 to 10	1
> 10	0

ISSUE 161

This revision updates language regarding the use of volatile chemicals in Type A1 BSCs in Sections 3.8.2.1, I-1.3.1.3, and I-1.11.8.2.1.

1. Class II Type A1 cabinets (formerly designated Type A): Cabinets that:

- maintain minimum average inflow velocity of 75 ft/min (0.38 m/s) through the work access opening;
- have HEPA/ULPA filtered downflow air that is a portion of the mixed downflow and inflow air from a common plenum (i.e., a plenum from which a portion of the air is exhausted from the cabinet and the remainder supplied to the total work area);
- may exhaust HEPA/ULPA filtered air back into the laboratory or to the environment through an external exhaust system connected to the cabinet with a canopy connection;
- have all biologically contaminated ducts and plenums under negative pressure or surrounded by negative pressure ducts and plenums; and
- **type A1 BSCs are not designed for work with volatile chemicals and are not suitable to be used for handling these chemicals.**

The highlighted language above was added and the language below was eliminated from the standard to avoid misapplication of the Type A1 cabinet. A1 cabinets are NOT to be used with volatile materials.

If using chemicals with toxic vapors, the unit shall be connected to an external exhaust system. Type A1 cabinets may be used for work with volatile chemicals if deemed appropriate by a chemical risk assessment

(refer to Section I-1.3.1.3).

NOTE — Type A1 BSCs manufactured prior to 2010 are not suitable for work with volatile chemicals due to the contaminated positive pressured plenums that are not surrounded by negative pressure plenums.

I-1.3.1.3, and I-1.11.8.2.1

Similar changes to those described for 3.8.2.1 were made to I-1.3.1.3 and I-1.11.8.2.1 to prohibit the use of Class II Type A1 cabinets with volatile materials.

ISSUE 162

Issue 162

This revision adds language regarding the feet per minute tolerance for the airflow velocity testing in Sections 6.9 and 6.17, and throughout Normative Annex 1.

9. Downflow and inflow velocities

1. The average downflow velocity (uniform downflow) or velocities (nonuniform downflow) and the calculated and measured average inflow velocities of the cabinet shall be set at the manufacturer's recommended nominal set points ± 2 ft/min (0.01 m/s) for testing or alternative set points as the test method requires. Subsequent production models of the test cabinets of the initial model and size conforming to Section 6.7 may also qualify when the inflow and average downflow velocity (or velocities, if so specified) operate within ± 5 ft/min (± 0.025 m/s) (see Section N-1.9) of the manufacturer's recommended nominal set points of the unit being tested.

*This section was modified from +/- 3 ft/min to +/- 2 ft/min. and minor rewording where indicated.
Previous wording ~~nominal set points ± 3 ft/min (0.015 m/s) for testing unless otherwise noted.~~*

ISSUE 162 (CONTINUED)

Issue 162

This revision adds language regarding the feet per minute tolerance for the airflow velocity testing in Sections 6.9 and 6.17, and throughout Normative Annex 1.

6.17 Air velocity stability

Air velocity stability shall be determined with the cabinet operating at the manufacturer's recommended nominal set points **± 2 ft/min (0.01 m/s)**

This section was modified from +/- 3 ft/min to +/- 2 ft/min.

Previous wording ~~Air velocity stability shall be determined with the cabinet operating at the nominal set point velocities ± 3 ft/min (0.015 m/s).~~

ISSUE 166

Issue 166

This revision updates language regarding the various uses of the term “dedicated” in Sections 3.8, I-1.3.1.4, and Table I-1.2.

I-1.3.1.4 Question four: If the BSC requires an exhaust system, is there an appropriate location for the cabinet and its ductwork?

— canopy-connected Types A and C1 require a consistent, low static pressure. While an exhaust duct **dedicated to that individual cabinet** is preferred, they may share a common exhaust system with other exhausted laboratory devices, if properly balanced;

*This section was modified from — canopy-connected Types A and C1 require a consistent, low static pressure. **While a dedicated exhaust system** is preferred, they may share a common exhaust system with other exhausted laboratory devices, if properly balanced;*

The new language is intended to specify that dedicated in this case implies “to that specific cabinet” as apposed to “to that building” or “to Class II BSCs. etc. This is only one example of where the standard further defines “dedicated” to mean “to that individual cabinet”

ISSUE 166 (CONTINUED)

Issue 166

This revision updates language regarding the various uses of the term “dedicated” in Sections 3.8, I-1.3.1.4, and Table I-1.2.

I-1.3.1.4 Question four: If the BSC requires an exhaust system, is there an appropriate location for the cabinet and its ductwork?

— Type B BSCs require a higher static pressure that must increase as their exhaust filters load. They must be on an exhaust duct and fan **dedicated to that individual cabinet**, and not be ganged with other Type B BSCs, or other exhausted laboratory devices requiring a lower static pressure (e.g., fume hoods, canopy-connected BSCs);

This section was modified from — Type B BSCs require a higher static pressure that must increase as their exhaust filters load. They must be on **a dedicated exhaust system**, and not be ganged with other Type B BSCs, or other exhausted laboratory devices requiring a lower static pressure (e.g., fume hoods, canopy-connected BSCs);

The new language is intended to specify that dedicated in this case implies “to that specific cabinet” as apposed to “to that building” or “to Class II BSCs. etc.

ISSUE 166 (CONTINUED)

Issue 166

This revision updates language regarding the various uses of the term “dedicated” in Sections 3.8, I-1.3.1.4, and Table I-1.2.

Table I-1.2 Characteristics of Type B1 and Type B2 BSCs

exhaust system type - Must have an exhaust duct and fan **dedicated to each BSC** - Must have an exhaust duct and fan **dedicated to each BSC**.

This section was modified from

*exhaust system type - Must have **dedicated ductwork and exhaust blower for each BSC** - Must have **dedicated ductwork and exhaust blower for each BSC**.*

The new language is intended to specify that dedicated in this case implies “to that specific cabinet” as apposed to “to that building” or “to Class II BSCs. etc. This is very moderate change to the language to be consistent with changes made elsewhere.

Note that similar changes are made elsewhere in this table.

ISSUE 167

Issue 167

This revision modifies language regarding the resistance to tipping, distortion and deflection of work surface under load tests in Section 6.8 and Normative Annex 1. It also removes former Figures 21 through 24.

8. Stability

The cabinet shall be designed and constructed to resist overturning and distortion under applied forces, resist deflection of the work surfaces under load, and resist tipping under workload.

1. Resistance to overturning

Cabinets shall conform to the requirements of UL 61010-117 or current edition, Section 7.3.

The sections above remain unchanged as posted here. The sections below were removed in their entirety to reduce the cost of cabinet accreditation. Stability testing is adequately covered in UL 61010-117.

2. Resistance to distortion

3. Resistance to deflection of work surface

4. Resistance to tipping

ISSUE 168

Issue 168

This revision adds language regarding ammonium carbonate in Section I-2.1.4.1.

I-2.1.4.1 Paraformaldehyde

CAUTION — All sources of hydrogen chloride must be removed from the cabinet before decontamination. Hydrogen chloride in the presence of formaldehyde, at ambient air conditions, will form the carcinogen Bis(chloromethyl)ether (BCME).

- a) Calculate the total volume of the cabinet by multiplying the height, width, and depth.
- b) Multiply the total volume of the cabinet by 0.30 g/ft³ (11 g/m³) of space to determine the gram weight of paraformaldehyde required [CHECK CONCENTRATION]. Determine the stoichiometric amount of **ammonium carbonate** [(NH₄)₂CO₃] or alternative to neutralize the resulting formaldehyde gas with ammonia gas. The ammonium carbonate should be weighed out so that it is 10% greater than the weight of paraformaldehyde used for the decontamination to ensure completion of the reaction.

*The words **ammonium carbonate** were added to any reference to the scientific notation [(NH₄)₂CO₃] throughout I-2.1.4.1. This is just one example.*

ISSUE 170

Issue 170

This revision updates language regarding the airflow smoke patterns test in Sections N-1.10.4.2 and N-5.4.4.2

N-1.10.4.2 **Sash** retention test

The visible aerosol or mist shall show smooth downward flow with no dead spots or reflux. No visible aerosol or mist shall escape from the cabinet.

N-5.4.4.2 **Sash** retention test

The visible aerosol or mist shall show smooth downward flow with no dead spots or reflux. No visible aerosol or mist shall escape from the cabinet.

*Both above changes were to rename this test “**Sash**” retention from “**View Screen**” retention for consistency.*

ISSUE 171

Issue 171

This revision updates the definitions of *Class II* and *decontamination* in Sections 3 and I-1-11. It also adds a footnote regarding HVAC in Sections 5.25.5 and I-1.3.1.

3.8.2 Class II: ...

All Class II cabinets are designed for work involving procedures assigned **to BSL-1, BSL-2, and BSL-3 facilities**. Class II BSCs may be used with procedures requiring BSL-4 facilities if used in a BSL-4 suit laboratory by a worker wearing a positive pressure protective suit.

*The **highlighted** section above was changed from **previous language below** for clarity. This same change was made in I-1.11*

All Class II cabinets are designed for work involving procedures assigned **to BSL 1-3 facilities**. Class II BSCs may be used with procedures requiring BSL-4 facilities if used in a BSL-4 suit laboratory by a worker wearing a positive pressure protective suit.

ISSUE 171(CONTINUED)

Issue 171

This revision updates the definitions of *Class II* and *decontamination* in Sections 3 and I-1-11. It also adds a footnote regarding HVAC in Sections 5.25.5 and I-1.3.1.

Section 5.25.5

²² **ASHRAE Laboratory Design Guide: Planning and Operation of Laboratory HVAC Systems, Second edition (2015) or latest edition. American Society of Heating, Refrigeration, and Air-Conditioning Engineers, Inc. 1791 Tullie Circle, NE, Atlanta, GA 30329. <www.ashrae.org>**

²³ Sheet Metal and Air Conditioning Contractors' National Association. 4201 Lafayette Center Drive, Chantilly, VA 20151. www.smacna.org

The **highlighted** footnote (22) was added. Footnote (23) is carried over.

ISSUE 171 (CONTINUED)

Issue 171

This revision updates the definitions of *Class II* and *decontamination* in Sections 3 and I-1-11. It also adds a footnote regarding HVAC in Sections 5.25.5 and I-1.3.1.

Section I-1.3.1

⁵⁷ Specifications taken from Military Specification MIL-F-51068 (Cancelled), Gasket Assembly.

[<www.defense.gov>](http://www.defense.gov)

The **highlighted** section above was changed from the **bold** section below

⁵⁷ Specifications taken from Military Specification MIL-F-51068 (Cancelled), Gasket Assembly.

[<www.defenselink.mil/pubs/>](http://www.defenselink.mil/pubs/)

NSF/ANSI 49: MARCH 2021 CHANGES

Changes Effective March 2021

The following is based on the changes outlined in the Foreword of NSF/ANSI 49

The topics were discussed in committee as issue numbers which form the basis of discussion.

- Issue 130, 141, 151, 153, 154, 155, 156, 157, 159, 160
- Only completed issues are discussed.

ISSUE 130

This revision updated language throughout the standard regarding use of the previously added terms total work area and usable work area.

- **N-1.6.1 Purpose**

- These tests determine whether aerosols will be contained within the cabinet, outside contaminants will not enter the cabinet usable work area, and aerosol contamination of other equipment in the cabinet will be minimized.

- **N-1.8.4 Acceptance**

- The average downward airflow velocity through the cross section of the unobstructed usable work area (with removable acceptable option components removed) at the level 4 inches (100 mm) above the bottom of the sash of cabinets meeting the requirements of Section N-1.6 shall be the values specified by the manufacturer.

- **I-1.3.1.1 Question one: What needs to be protected?**

- ...The Class I BSC will protect the operator and the lab, however, because room air constantly washes over the usable work area, the product is exposed to airborne contaminants.

ISSUE 130 (CONTINUED)

I-1.7.1.5 Effective usable work area layouts will minimize reach to avoid neck and shoulder stress and fatigue. Rotating tables are available to minimize reach.

I-1.7.4 BSC techniques

- ...and (2) airborne contamination generated in the usable work area is controlled by flow of airstreams in a top-to-bottom direction;

I-1.7.5 BSC start up procedure

- ...This will minimize the shedding of skin flora into the usable work area and also protect hands and arms from viable microbial contamination.

I-1.11.2.3 downflow velocity profile:

- ...Airflow velocities and the average of the airflow through the usable work area may be calculated as a whole (uniform) or may be separated into two or more adjoining areas (zoned) with averages calculated for each zone.

Note there are 92 instances of the use of “total work area”. Too numerous to list here.

ISSUE 141

This revision updated multiple definitions in Section 3

- Some definitions moved from other sections into the definitions section 3.
- **3.8 cabinet classification:** Although this standard covers only Class II BSCs, Class I and Class III cabinets are currently defined and known to be commercially available.
 - This definition was modified.
- Other minor changes throughout the standard

ISSUE 151

This issue modified language in Section N-1.6 regarding the occurrence of positive control plates.

ISSUE 153

This revision updated language in Sections 3.19.1, 6.11, 6.15, and N-1.6.5.1 regarding use of the term *smoke*.

- The word “smoke” was removed from the standard because the visible medium used is not truly a smoke.
 - 3.19.1 Changed smoke particles to aerosol particles
 - 6.11, 6.15 Changed “Airflow smoke patterns” to “Airflow pattern test” and removed any reference to “smoke”
 - 6.1.6.5.1 Changed reference to “smoke” to “visible medium”

ISSUE 154

This revision modified language in Sections N-1.8, N-1.9, N-5.2, N-5.3, and N-5.12 regarding air velocity measurements and the field certification label.

- N-1.8 added “The anemometer probe shall not be hand held for any downflow velocity measurements.”
- N-1.9 added “The anemometer probe shall not be hand held”
- N-5.2 added “The anemometer probe shall not be hand held for any downflow velocity measurements.”
- N-5.3 added “The anemometer probe shall not be hand held. Acceptable methods include a ring-stand and clamp, manufacturer supplied probe holder or when the BSC manufacturer has made provisions for accurately locating the anemometer sensing element by resting the probe on the exhaust collar and a second lip for stability.
- N-5.12 added “Technician NSF listing number, if applicable.” to the Field certification label requirements.

ISSUE 155

This revision updated language in Section 5.25.6 regarding the warm up period.

- This sentence was changed as follows: When starting the cabinet blowers from a dead stop, the inflow alarm must activate a visual indication until the cabinet either enters into a visually indicated warm up period ~~not to exceed 2 minutes~~ or the appropriate inflow velocity is achieved to ensure proper alarm system function.

ISSUE 156

This revision affirms language in Section N-1.6 regarding the placement of biological test control plates.

- Added to N-1.6.5.1.4.b) “A visible aerosol or mist test may be performed to determine where the test organism will be best captured in the areas described above, allowing unnecessary control plates to be eliminated.”

ISSUE 157

This revision further clarifies language in Section N-1.6 regarding biological tests.

ISSUE 159

This revision modifies language regarding electrical safety in Section 6.14.

Added “For the purposes of this requirement, an “authorized testing laboratory” shall be either a US Occupational Safety and Health Administration (OSHA) Nationally Recognized Testing Laboratory (NRTL) or an International Electrotechnical Commission for Electrical Equipment CB Testing Laboratory (CBTL).”

ISSUE 160

3.8.2

- This sentence was eliminated “Class II BSCs provide the microbe-free work environment necessary for cell culture propagation and also may be used for the formulation of nonvolatile antineoplastic or chemotherapeutic drugs.”

ISSUE 160 CONTINUED

This revision adds clarifying language to definitions in Sections 3.5.3 and 3.8.2.

- **3.8.3 Class III:** The Class III BSC was designed for work with highly infectious microbiological agents and other hazardous operations. It provides maximum protection for the environment and the worker. ~~It is a gas-tight (no leak greater than 1×10^{-7} mL/s with 1% test gas at 3 inches (750 Pa) pressure water gauge)21 enclosure with a viewing window that is secured with locks, or requires the use of tools to open, or both.~~ Access for passage of materials into the cabinet may be through any of the following: a dunk tank that is accessible through the cabinet floor, a double-door pass-through box that can be decontaminated between uses, integrated double door autoclaves, and portable docking stations with double sealing connecting mechanisms that can be decontaminated between uses. Reversing that process allows materials to be removed from the Class III BSC. Both supply and exhaust air are HEPA/ULPA filtered. Exhaust air must pass through two HEPA/ULPA filters in series, before discharge to the outdoors. Airflow is maintained by an exhaust system exterior to the cabinet, which keeps the cabinet under negative pressure according to manufacturer design pressure criteria. ~~Sometimes because of laboratory conditions an optional exhaust fan may be required. This exhaust fan should generally be kept separate from the exhaust fans of the facility ventilation system. If a cabinet exhaust system is required it should be equipped with an appropriate alarm system which both notifies the cabinet user and shuts down the cabinet exhaust system in the event of a facility exhaust system failure.~~ This exhaust fan should generally be kept separate from the exhaust fans of the facility ventilation system and be equipped with an appropriate alarm system which both notifies the cabinet user and shuts down the cabinet exhaust system in the event of a facility exhaust system failure.

ISSUE 160 CONTINUED

3.8.3 Class III:

This paragraph was modified as noted:

- Rubber gloves / sleeves or equivalent glove material, are sealed to ports in the cabinet and allow direct manipulation of the materials isolated inside. The glove material shall be compatible for use with the materials being used in the cabinet. The exhaust system for the cabinet shall provide inflow to the cabinet arm port in the event of a rubber glove / sleeve breach. ~~The minimum breach velocity shall be measured with a hot wire in the middle of the arm port and shall be no less than 100 ft/min (0.51 m/s). It is not a requirement for the work area to be free of turbulence or cross contamination.~~

NSF/ANSI 49: MARCH 2020 CHANGES

Changes Effective March 2020

The following is based on the changes outlined in the Foreword of NSF/ANSI 49

The topics were discussed in committee as issue numbers which form the basis of discussion.

- Issue 54, 82, 92, 120, 122, 127, 130, 133, 136, 138, 139, 140, 146, 147, 148, 149
- Only completed issues are discussed.

ISSUE 54

This revision affirms new and updated language in Annex N-5 (formerly Annex F) concerning the use of the Secondary method for measuring airflow.

ISSUE 82

This revision affirms new language in Annex I-1 (formerly Annex E) regarding the term *percent recirculation*.

ISSUE 92

This revision affirms new and updated language in Section 3 and Annex N-5 (formerly Annex F) regarding canopy field testing.

ISSUE 120

This revision affirms new language in Section 3 regarding the addition of the newly proposed term *plenum*.

ISSUE 122

This revision affirms new and revised language in Annex N-5 (formerly Annex F) regarding the Certification Label.

ISSUE 127

This revision affirms revised language regarding the use of the term *NOTE*.

ISSUE 130

This revision affirms revised language regarding the definition of the term *work area*.

ISSUE 133

This revision affirms revised language in Section 5 regarding the data plate.

ISSUE 136

This revision affirms revised language in Annex N-1 (formerly Annex A) and Annex N-5 (formerly Annex F) regarding the sash seal smoke test.

ISSUE 138

This revision affirms revised language regarding the range of measurement for vibration frequency.

ISSUE 139

This revision affirms revised language in Annex I-1 (formerly Annex E).

ISSUE 140

This revision affirms revised language in Annex N-1 (formerly Annex A) regarding accuracy requirements for the manometer used for the pressure decay and motor blower performance tests.

ISSUE 146

This revision addresses inconsistencies of incubation times and temperatures during the various biological tests in Annex N-1.

ISSUE 147

This revision affirms revised language in Annex N-1 (formerly Annex A) regarding filter porosity for filtering impinger water.

ISSUE 148

This revision affirms revised language in Annex N-1 (formerly Annex A) regarding the confirmation requirements for the Cross Center test.

ISSUE 149

This revision affirms revised language in Section 2 regarding Normative References.

NSF/ANSI 49: MARCH 2020 CHANGES

Annexes

Annexes	
Previously known as:	Now known as:
Annex A	Normative Annex 1 (N-1)
Annex B	Normative Annex 2 (N-2)
Annex C	Normative Annex 3 (N-3)
Annex D	Normative Annex 4 (N-4)
Annex E	Informative Annex 1 (I-1)
Annex F	Normative Annex 5 (N-5)
Annex G	Informative Annex 2 (I-2)
Annex H	Informative Annex 3 (I-3)
Annex I	Informative Annex 4 (I-4)
Annex J	Informative Annex 5 (I-5)
Annex K	Informative Annex 6 (I-6)

NSF/ANSI 49: JANUARY 2019 CHANGES

The following is based on the changes outlined in the Foreword of NSF/ANSI 49

The topics were discussed in committee as issue numbers which form the basis of discussion.

- Issue 47, 59, 77, 105, 108, 109, 110, 111, 112, 115, 117, 118, 121, and 125
- Only completed issues are discussed.

ISSUE 47

Revised and added new language in Annex A regarding the cross contamination test procedure.

ISSUE 59

Revised and added new language for the airflow alarm requirement language for all cabinet types.

ISSUE 77

Added language to Annex A regarding confirmation tests when there is a test failure.

ISSUE 105

Added new language in Annex E regarding risk assessment of biosafety cabinets exhaust system pressurization in the event of an exhaust system failure.

ISSUE 108

Added new language in Annex F regarding exhaust airflow alarms

ISSUE 109

Revised and added new language referring to the term Certification throughout this standard.

ISSUE 110

Revised language in Section 1.3 regarding the identification of major software changes to biosafety cabinets.

ISSUE 111

Added new and revised language to Section 3 regarding the terms visible and viewable, and added provisions for the optional use of digital data plates to Section 5.

ISSUE 112

Revised and added new language in Annex G regarding the generation and dispersion of decontamination gas.

ISSUE 115

Added new language in Section 6.10 regarding inflow requirements for Type C1 Cabinets.

ISSUE 117

Updates the language throughout this Standard regarding the use of the term “and/or”.

ISSUE 118

Revised language in Section 5.25.3 and Annex F regarding the exhaust alarm in Type B biosafety cabinets.

ISSUE 121

Revised language in Section 5.26.2 regarding electrical wiring.

ISSUE 125

Revised language in Section 5.26.2 regarding the term “running power”.

NSF/ANSI 49: MARCH 2017 CHANGES

Changes Effective March 2017

Note: That in addition to this summary, CETA published an applications guide to provide a more detailed explanation including rationale for the changes.

The following is based on the changes outlined in the Foreword of NSF/ANSI 49

The topics were discussed in committee as issue numbers which form the basis of discussion.

- Issue 45, 56, 73, 76, 78, 79, 81, 86, 88, 90, 96, and 99
- Only completed issues are discussed.

ISSUE 45

Changes were made to Annex G covering the addition of Vaporized Hydrogen Peroxide as a decontamination agent in Biosafety Cabinetry, as well as clarifying the use of Chlorine Dioxide.

ISSUE 56

Language in section 6.14 was added regarding certification to IEC 61010-1 or a national standard based on it.

ISSUE 73

Updated/unified all cabinet definitions, Sections 5, 6, Annex A and F for the type C1.

ISSUE 76

Updated language was added referencing audible and visual alarms.

ISSUE 78

Metric conversions throughout the Standard were reviewed and updated.

ISSUE 79

Decontamination procedures in Annex G were updated.

ISSUE 81

Updated Annex E. Included all definitions from section 3, updated risk assessment, cabinet placement, cleaning procedures, usage procedures and added type C1.

ISSUE 86

Updates were made in Annex A to the incubation time and temperature prior to the micro check.

ISSUE 88

Section 5.32 was updated regarding cabinet height and width.

ISSUE 90

Redundancies about sliding sash alarms in subsections 5.19.4 and 5.25.1 were removed.

ISSUE 96

Language involving preparation of the spore suspension of *Bacillus atrophaeus* 9372 in Annex A was updated.

ISSUE 99

Update Soap bubble leak test procedure in Annex A to require manufacturers to provide metal seal plates and limits NSF test personnel time to seal leaks found.

NSF/ANSI STANDARD 49: FEBRUARY 2015 CHANGES

Changes effective February 2015

Note that in addition to this summary, CETA published an applications guide to provide a more detailed explanation including rationale for the changes.

- Issues 48, 49, 50, 51, 52, 53, 55, 60, 61, 72

ISSUE 48: MOTOR STABILITY TEST PROCEDURE

This revision adds a motor stability test procedure for motor speed control systems.

ISSUE 49: SEALANT USE LANGUAGE

This revision updates the sealant use language in Annex H, Section H.6.

- Updated recommended sealants for the BSCs.
 - SAE AMS-S-8802 or equivalent
 - Dow Corning RTV 732, 781, 734, or RRTV 112 self leveling or equivalent

ISSUE 50: FANS

This revision affirms new language regarding the type of fans used in biosafety cabinets.

ISSUE 51: BIOSAFETY CABINET BLOWER STARTUP

This revision affirms new language regarding the type of biosafety cabinet blower startup.

ISSUE 52: DOP PORT LOCATION

This revision clarifies details surrounding the DOP port location in section 5.22

- Moved the challenge port from the top of the BSC on the outside to the inside under the work surface. The intent is to ensure that the challenge port is opened in negative pressure.

ISSUE 53: DEFINITIONS

This revision adds definitions to clarify biosafety cabinet shell penetrations and cable ports with consideration given to service technicians and cabinet users relating to safety.

ISSUE 55: INSTRUMENTATION LANGUAGE

This revision updates the instrumentation language.

ISSUE 60: AIRFLOW GRID LANGUAGE

This revision updates the airflow grid language in sections A.8.3.1 and A.8.3.2, and the related figure A15

- Clarifies minimum spacing and number of readings determination.
- Clarifies different zone dimensions for cabinets with sloped sash for different sash height

ISSUE 61: SASH POSITION

This revision updates the language in sections 5.19.4 and 5.25.1 to include a section requiring the use of a sash position “too low” alarm.

- Now must alarm for more than 1” above and below

ISSUE 72: FIGURES

This revision updates multiple figures throughout the Standard to improve clarity.

NSF/ANSI STANDARD 49: JULY 2012 CHANGES

Changes effective July 2012

- Issues 44, 46

ISSUE 44: CABINETS LESS THAN 3 FEET

Updates the language in the Standard to include a test method for biosafety cabinets with an interior sidewall dimension of three feet or less.

ISSUE 46: DIRECT INFLOW MEASUREMENT (DIM) & EXHAUST CLEARANCE

Updates the Standard to include a reference to non-back pressure compensated readings used in a Direct Inflow Measurement (DIM) in Annex A, Annex B, and Annex F

Updates the language in Annex A, Annex E and Annex F for the 12 in (30 CM) clearance requirement used for measuring an exhaust HEPA filter

NSF/ANSI 49: NOVEMBER 2010 CHANGES

Changes effective November 29, 2010

The following summary of changes is based on the changes outlined in the Foreword

The issues were discussed in committee as issue numbers which form the basis of discussion.

- Issue 23, 29
- Only completed issues are discussed.
- Note that these two simple changes are addendum to the other 2010 changes

ISSUE 23: B2 DEFINITION

Errors in the September 2010 version were corrected. B2 Cabinet description is now correct.

ISSUE 29: UNIFORM AND ZONED DOWNFLOW

1. Uniform downflow cabinets

- No change to this section

2. Non-uniform downflow (zoned) cabinets

- The grid must have equidistant spacing
- Each zone must have at least 7 points within it
- The distance between test points shall not be less than 4" nor more than 8" inches apart
- The area defined by the perimeter of the test points must equal at least 30% of the total area of the plane in which the readings are taken.
- Each zone shall be taken at least 6 in away from the walls and sash enclosing the work area.

NSF/ANSI 49: SEPTEMBER 2010 CHANGES

Changes effective September 2010

The following summary of changes is based on the changes outlined in the Foreword

The issues were discussed in committee as issue numbers which form the basis of discussion.

- Issue 23, 24, 37/38, 41
- Only completed issues are discussed.

One additional item was added and not listed as with an issue number.

- Annex K

ISSUE 23: HARD DUCTING CABINETS

2. “The canopy connection type of BSC exhaust connection is required for externally vented Class II, Type A1 or A2 BSCs.”

- Minor changes to wording from previous version of standard.

F.7.3.3 “Direct connected A1 or A2 BSCs shall not be considered in compliance with the Standard.”

- New language intended to end “grace period” for conversion of old installation direct connections for type A1 and A2 cabinets.
 - Strong recommendation since 2002 now becomes a requirement. All externally vented A cabinets now must be connected with a canopy exhaust connection or they cannot be certified to be in compliance with the standard.

ISSUE 23: HARD DUCTING CABINETS

F.7.3.3 “Using a visible medium source positioned to demonstrate containment of BSC exhaust by the canopy, reduce the external exhaust until the alarm signals audibly. The alarm shall sound before visible canopy containment is lost.”

- New language for testing canopy connections.
- It should be noted that new language for this is anticipated that will be consistent with other alarm requirements (within 15 seconds).

ISSUE 23: ULPA FILTERS

HEPA references now are HEPA/ULPA

ISSUE 24: ALARMS

5.23.4 “Any Type A1 or A2 cabinet when canopy connected shall have audible and visual alarm indication notifying the user of a potential loss in canopy containment”

ISSUE 37/38: ILLUSTRATIONS

New illustrations throughout the document

ISSUE 41: IEC 61010

“The standard was revised to be more inclusive of markets outside North America by modifying section 6 – Performance of the standard.”

- **6.14 Electrical safety**

“The cabinet shall be tested by a National Recognized testing Laboratory (NRTL) for compliance to the requirements of the current edition of any national standard that is based on IEC 61010-1. Compliance is demonstrated by cabinet listing, i.e. UL, CSA, or IECEE CB Scheme certificate.”

ANNEX G: AMMONIUM CARBONATE (NOT A NUMBERED ISSUE)

1. “The ammonium carbonate should be weighed out so that it is 10% greater than the weight of paraformaldehyde used for the decontamination to ensure completion of the reaction.”

- Changed from Ammonium bicarbonate to reflect Dr. Luftman’s work published in “Applied Biosafety 10(2)-2005.”
 - Bicarbonate should be 1.6:1 according to this article.

ANNEX K

A “Protocol for the Validation of Alternative Biosafety Cabinet Decontaminating Methods and Agents” was added in the form a new annex (Annex K).

- This is an informative annex

NSF/ANSI 49: JUNE 2009 CHANGES

Changes effective June 2009

The following summary of changes is based on the changes outlined in the Foreword

The issues were discussed in committee as issue numbers which form the basis of discussion.

- Annex G, Issue 15, 28, 30, 34, 35, 36
- Only completed issues are discussed.

ANNEX G

New section in 2008 largely unchanged.

ISSUE 15: AEROSOL INTRODUCTION POINT

A2.3.1a) “The manufacturer shall determine the aerosol introduction point that provides the most uniform distribution”

- Same as 2008

New language for 2009

- The **location of the aerosol introduction point** shall be clearly described or indicated in a manner readily available to the certifier.
 - On the cabinet data plate (or)
 - With the electrical schematic if the schematic is affixed to the cabinet

ISSUE 28: ILLUSTRATIONS

Illustrations throughout the document were updated to take advantage of modern technology.

- Thanks to Jim Hunter, Labconco Corporation

ISSUE 30: BIOLOGICAL VS. BIOSAFETY

The term “biosafety” is now used throughout the standard. The title was changed previously but there were some cases where old terminology remained in the 2008 version.

ISSUE 34: ACCEPTANCE STATEMENTS

Correct the acceptance statements in Annex F for consistency with Annex A.

- F2.4 Format change. Use of word “shall”
- F3.4 Use of word “shall”
- F4.4 No change
- F5.4 No change
- F6.1.5 Use of word “shall”
- F7 N/A
- F8 N/A
- F9.4 Use of word “shall”
- F10.4 Use of word “shall”
- F11.4 Use of word “shall”

ISSUE 35: MULTIPLE ISSUES

3.26 Added a definition for w.g. (water gauge)

- “Another common name for inch of water column. The word “gauge” after a pressure reading indicates that the pressure stated is actually the difference between the absolute or total pressure and the air pressure at the time of the reading”.

UL References were updated throughout.

1. First reference to NSF/ANSI 49-2002 changed to NSF/ANSI 49 to clarify.

- With exception of downflow velocity test, **all** cabinets will be field tested to **current** version of Annex F regardless of manufacture date.

ISSUE 36: PLENUM DESIGN

- 5.4 All biologically contaminated ducts and plenums in Type A1, A2, B1, and B2 cabinets shall be maintained under negative pressure or enclosed within a negative pressure zone.**
- Match the definition for a type A1 cabinet as in the 2008 version of NSF/ANSI 49. Removed contradictory language.
 - A1 cabinets can no longer have positive pressure contaminated plenum.

NSF/ANSI STANDARD 49: APRIL 2008 CHANGES

Changes effective April 2008

- Issue 12, 13, 14, 16, 17, 18, 19, 20, 22, 25, 26, 31

ISSUE 12: HEPA FILTERS

3.14* High efficiency air filters (for use in class II biosafety cabinets)

- Added **(for use in class II biosafety cabinets)** to the title to make it clear that the language in the standard only addresses filters as they pertain to BSCs.

*Note that this is NOT section 3.13 as stated in the foreword.

ISSUE 13, 21: CABINET DESCRIPTIONS - TYPE A

Change for Class II Type A1

- All biologically contaminated plenums under negative pressure or surrounded by negative pressure ducts and plenums.
 - Positive pressure exterior plenums are no longer allowed on any Class II BSC.

1. Tests directly related to containment – cabinet integrity test

- Old A1 cabinets only
- New positive pressure plenum cabinets only

ISSUE 14: CONCURRENT BALANCE VALUE

Added definition for CBV (3.10)

- For all direct-connected BSCs
- Compares the primary (DIM) value to a duct traverse.
 - Traverse per ASHRAE std. 111-2008
- Requirements for listing now specified
 - Exhaust volume @ filter load value (pressure)
 - Allowance for filter loading added to measured value
 - B1 0.3" w.c.
 - B2 0.7" w.c.

ISSUE 16: STANDARD TITLE

Biosafety Cabinetry: Design, Construction, Performance, and Field Certification

- Remove reference to Class II to include all types of cabinets

ISSUE 17, 20: INTERLOCKS FOR TYPE B CABINETS

Section F.7.3.2

- **Changed requirement for interlock verification to Type B cabinets from Type B2 cabinets.**
 - Significant impact for B1 cabinet testing!

ISSUE 18: DOWNFLOW VELOCITY MEASUREMENT

1. clarification language

- Downflow velocity readings shall be taken 4" above the bottom edge of the window only when so stated on the manufacturers data plate or when the manufacturers data plate label indicates the cabinet was listed to NSF 49-2002 or later.

ISSUE 19: SOUND LEVEL MEASUREMENTS - ANNEX F

F.11.4 Acceptance (sound level)

- Cabinet passes when overall noise level does not exceed 70 dbA when ambient is not greater than 60 dbA.
 - When ambient exceeds 60 dbA, correction curves are used.

Difference between total and background sound readings in dbA	Number to subtract from total to yield corrected noise level
0-2	Reduce background levels
3	3
4-5	2
6-10	1
>10	0

ISSUE 26: REPORTED VALUES

Certification Report - F.12.2

- A certification report that will carry the language “**certified in accordance with the NSF annex F**” or any similar language shall, at a minimum, include the following:
 - BSC Model Number
 - BSC Serial Number
 - BSC Location
 - **BSC Venting Information**
 - Ducted or not ducted
 - Type of connection (canopy, direct, or none)
 - Type of BSC
 - Test equipment used for each test
 - Manufacturer, model, serial number, calibration date
 - **Specific test data as detailed in annex F**
 - **Acceptance criteria for each test**
 - **Printed** name of certification technician
 - Retest date

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Downflow Velocity for Uniform Downflow Cabinets – F.2.3.1**
 - Individual velocity readings in the applicable grid
 - Overall average of velocity readings
 - Minimum velocity reading
 - Maximum velocity reading
 - Acceptance criteria for average airflow velocity
 - Acceptance criteria for airflow velocity uniformity
 - Name of test (Uniform Downflow Velocity Test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Downflow velocity for Non-Uniform Downflow Cabinets** *for each zone* – **F.2.3.2**
 - Individual velocity readings in the applicable grid
 - Overall average of velocity readings
 - Minimum velocity reading
 - Maximum velocity reading
 - Acceptance criteria for average airflow velocity
 - Acceptance criteria for airflow velocity uniformity
 - Name of test (Uniform Downflow Velocity Test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Intake Velocity using the Direct Measurement Method – F.3.3.2**
 - Individual volume readings
 - Overall average of the volume
 - Calculated inflow volume
 - Work access opening area
 - View screen opening height
 - Correction factor used (if applicable)
 - Acceptance criteria for average inflow volume
 - Acceptance criteria for calculated inflow velocity
 - Inflow velocity test method
 - Name of test (Inflow velocity test)

ISSUE 26: REPORTED VALUES

Reported Values

Annex F: Requirements for “Reported Values” are now delineated for each test in Annex F.

- **Intake Velocity using a thermal anemometer to measure exhaust velocity to determine inflow velocity – F.3.3.3.1**
 - Individual exhaust velocity readings
 - Overall average of the exhaust velocity readings
 - Calculated exhaust volume
 - Calculated inflow volume
 - Exhaust opening dimensions
 - Exhaust opening effective area
 - Work access opening area and dimensions
 - View screen opening height
 - Correction factor used (if applicable)
 - Acceptance criteria for calculated inflow velocity
 - Acceptance criteria for calculated inflow velocity
 - Inflow velocity test method
 - Name of test (Inflow velocity test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Intake Velocity using a thermal anemometer to measure velocity through a constricted access opening to determine average inflow velocity – F.3.3.3.2**
 - Individual constricted velocity readings
 - Overall average of the inflow velocity readings
 - Calculated inflow volume
 - Work access opening dimensions and area
 - Correction factor used (if applicable)
 - Acceptance criteria for average inflow velocity
 - Inflow velocity test method
 - Name of test (Inflow velocity test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Intake Velocity using a thermal anemometer to measure velocity through the access opening to determine average inflow velocity (B1) – F.3.3.3.3**
 - Individual velocity readings
 - Overall average of the inflow velocity readings
 - Calculated inflow volume
 - Work access opening dimensions and area
 - Correction factor used (if applicable)
 - Acceptance criteria for average inflow velocity
 - Inflow velocity test method
 - Name of test (Inflow velocity test)

ISSUE 26: REPORTED VALUES

Reported Values

Annex F: Requirements for “Reported Values” are now delineated for each test in Annex F.

– Intake Velocity using an anemometer and pitot tube – F.3.3.3.4

- Individual duct velocity readings
- Overall average of the duct velocity readings
- Calculated exhaust volume
- Duct size shape and area
- Calculated inflow volume
- Work access opening dimensions and area
- Dimensions and area of the supply velocity measurement locations used to determine supply volume
- Individual supply velocity readings (not to be confused with downflow velocities)
- Calculated supply velocity and volume
- Calculated inflow velocity and method used for calculations
- Correction factor used (if applicable)
- Acceptance criteria for calculated inflow velocity
- Inflow velocity test method
- Name of test (Inflow velocity test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Airflow Smoke Pattern Test – F.4.3**
 - Name of each test
 - Pass or fail for each test

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **HEPA Filter Leak Test – F.5.3.1**
 - Upstream Aerosol Challenge Concentration
 - Method used to report concentration (measured or calculated)
 - Maximum leak penetration in percent
 - Method used (scanned or Probe tested)
 - Name of test (HEPA filter leak test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Pressure decay / soap bubble test – F.6.1.3**
 - Pressure range maintained during test
 - Pass or fail
 - Name of test (Pressure Decay Test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Site installation assessment tests – F.7**
 - Reported values are not specified in standard

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Lighting Intensity Test – F.9.3**
 - Individual background readings
 - Individual lighting intensity readings
 - Average background intensity
 - Average lighting intensity
 - Acceptance criteria
 - Pass or fail
 - Name of test (Lighting intensity test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Vibration Test – F.10.3**
 - Unit “On” vibration reading
 - Background vibration reading
 - Net vibration
 - Pass or fail
 - Name of test (Vibration Test)

ISSUE 26: REPORTED VALUES

Reported Values

- Annex F
 - Requirements for “Reported Values” are now delineated for each test in Annex F.
 - **Noise Level Tests – F.11.3**
 - Unit “On” sound level reading
 - Background sound level reading
 - Net sound level
 - Pass or fail
 - Name of test (Noise level test)

ISSUE 22, 25: ANNEX G

Informative annex expanded from formaldehyde decontamination to include:

- G.1 Biosafety Consultation prior to BSC purchase
- G.2 Risk Assessment procedures
- G.3 Cabinet selection
- G.4 Prior to the Purchase
- G.5 Inspection (of new cabinets)
- G.6 Moving a BSC
- G.7 Decontamination procedures
 - Expanded to include Chlorine Dioxide procedures
- G.8 HEPA Filter Disposal Procedures
- G.9 Lifespan of BSCs
- G.10 Decommissioning process

ISSUE 31: HELIUM AND SULFUR HEXAFLUORIDE TESTS

Annex J created as informational annex because there are no longer any exterior positive pressure contaminated plenum designs. Material removed from Annex A.

- J.1 Helium Leak Test
 - Previously section A.1
- J.2 Sulfur Hexafluoride (SF₆) Leak Test
 - Previously section A.2

SUMMARY

Both minor and major changes were made to NSF/ANSI 49-2008. The most substantive changes were elimination of Type A cabinets with exterior positive pressure contaminated plenums, listing requirements for the Concurrent Balance Value, inclusion of chlorine dioxide decontamination procedures, and requirements for minimum reported values for field certification. Minor changes were made in 2009. The most substantive change made in 2010 is strengthening of language relating to external venting of class II Type A cabinets and specific new guidance for testing canopy connections.